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Dental Implant Abutment Geo Customized Abutment

GM-IFU-KR-CA, 2026.10.29 (Rev.0)

Intended use
The product is an implant abutment that connects to a dental implant to support prostheses, such as artificial teeth, used to restore patients' masticatory function.

Usage
A. Before use
1. Check if the product is normal and free from any foreign substance before use. Do not use the product if the expiration date has passed.
2. Use the standardized surgical instruments(driver, wrench) that properly fit the product, and do not apply excessive force.
3. The product should be used only by a skilled clinician or authorized person. It is intended for single use only and must not be reused under any circumstances.
4. All abutments and screws must be sterilized at 121℃ for 30 minutes prior to use.
5. Before placing the superstructure into the patient's oral cavity, the clinician must confirm the osseointegration status of the implanted fixture through radiographic evaluation and percussion tests, and proceed only if appropriate.
6. The clinician should be fully familiar with the surgical procedures, clinical indications, and precautions associated with this product and its instruments.
7. The patient's oral condition should be carefully evaluated to ensure that there are no biological or physiological factors that may adversely affect the surgical outcome.
8. This product is a customized device and must be appropriately designed and manufactured, taking into full consideration the patient's oral condition, implant position, and occlusal relationship.

B. How to use
1. Design a customized abutment considering the patient's oral condition, implant position, and occlusion.
2. Acquire the patient's oral information using a digital scanner or impression taking, and perform the design using CAD software.
3. Generate machining data using CAM software, and fabricate the abutment using a CAD/CAM milling machine.
4. After removing the healing abutment, connect the appropriate customized abutment to the fixture using the recommended tightening torque.
5. The abutment is designed and manufactured in various connection structures and configurations

C. How to store and manage after use
This product is intended for **single use only.** It should **not be reused.**

Precautions for Use
A. Precautions for Use
1. This product must not be used for any purpose other than its intended use.
2. Before use in the patient's oral cavity, the clinician must confirm the osseointegration status of the fixture through radiographic and clinical evaluation.
3. As this is a non-sterile product, it must be sterilized at 121℃ for 30 minutes prior to intraoral use.
4. Use standardized surgical instruments (driver, wrench) and do not apply excessive force.
5. The product should be used only by a skilled clinician or authorized person and is intended for single use only. Do not reuse under any circumstances.
6. Check the expiration date and packaging condition before use. Do not use the product if the packaging is damaged or the expiration date has passed.
7. The abutment must not be designed or machined beyond the approved design range. Excessive machining or structural modifications may adversely affect the performance of the product.
8. This product is a customized device. Errors in the design or manufacturing process may lead to poor fit, loosening, or fracture. Therefore, sufficient verification must be performed at each stage.

B. Side effects
Local allergic reactions may occur.

C. Warning
1. The service life of the prosthetic components is difficult to predict and is considered semi-permanent. Due to various biological, mechanical, and physicochemical factors, the durability of the implanted components cannot be guaranteed. Patients must be fully informed of this.
2. Inadequate initial fixation, potential infection, trauma, or the transmission of loads to the implant may result in mechanical loosening.
3. Improper fixation may lead to implant damage or fracture.
4. This product is a customized device and must not be used for any patient other than the patient for whom it was designed.

Storage
Store at room temperature.